

Heat-stable enterotoxin receptor

Catalog No: #AP84530



Package Size: #AP84530-1 50ug #AP84530-2 100ug #AP84530-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat-stable enterotoxin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	STA receptor, Guanylyl cyclase C, GC-C, Intestinal guanylate cyclase, Guc2c
Accession No.	Uniprot ID:Q3UWA6
Uniprot	Q3UWA6
GeneID	14917;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gucy2c

Note: This product is for in vitro research use only