

# Prostatic glandular kallikrein-6

Catalog No: #AP84533



Package Size: #AP84533-1 50ug #AP84533-2 100ug #AP84533-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Prostatic glandular kallikrein-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Glandular kallikrein-8,rGK-8,P1 kallikrein,Tissue kallikrein,Klk-8,Klk8
Accession No.	Uniprot ID:P36374
Uniprot	P36374
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Klk6

Note: This product is for in vitro research use only