

Ladinin-1

Catalog No: #AP84536



Package Size: #AP84536-1 50ug #AP84536-2 100ug #AP84536-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ladinin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lad-1,Linear IgA disease antigen,LADA,LAD
Accession No.	Uniprot ID:O00515
Uniprot	O00515
GeneID	3898;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LAD1

Note: This product is for in vitro research use only