

Type II iodothyronine deiodinase

Catalog No: #AP84537



Package Size: #AP84537-1 50ug #AP84537-2 100ug #AP84537-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Type II iodothyronine deiodinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	5DII,DIOII, Type 2 DI, Type-II 5'-deiodinase
Accession No.	Uniprot ID:Q9Z1Y9
Uniprot	Q9Z1Y9
GeneID	13371;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Dio2

Note: This product is for in vitro research use only