

Thyroxine 5-deiodinase

Catalog No: #AP84538



Package Size: #AP84538-1 50ug #AP84538-2 100ug #AP84538-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thyroxine 5-deiodinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	1.21.99.3,5DIII,DIOIII,Type 3 DI,Type III iodothyronine deiodinase
Accession No.	Uniprot ID:Q6QN11
Uniprot	Q6QN11
GeneID	414378;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DIO3

Note: This product is for in vitro research use only