

## Heat shock protein beta-1

Catalog No: #AP84539



Package Size: #AP84539-1 50ug #AP84539-2 100ug #AP84539-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Heat shock protein beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HspB1,Growth-related 25 kDa protein,Heat shock 25 kDa protein,HSP 25,Heat shock 27 kDa protein,HSP 27,p25,Hsp25,Hsp27
Accession No.	Uniprot ID:P14602
Uniprot	P14602
GeneID	15507;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Hspb1

Note: This product is for in vitro research use only