

Interleukin-31 receptor subunit alpha

Catalog No: #AP84542



Package Size: #AP84542-1 50ug #AP84542-2 100ug #AP84542-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-31 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-31 receptor subunit alpha, Cytokine receptor-like 3, GLM-R, hGLM-R, Gp130-like monocyte receptor, Gp130-like receptor, ZcytoR17, CRL3, GPL, UNQ6368/PRO21073/PRO21384
Accession No.	Uniprot ID:Q8NI17
Uniprot	Q8NI17
GeneID	133396;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL31RA

Note: This product is for in vitro research use only