

Inhibin beta B chain

Catalog No: #AP84545



Package Size: #AP84545-1 50ug #AP84545-2 100ug #AP84545-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inhibin beta B chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Activin beta-B chain
Accession No.	Uniprot ID:P04088
Uniprot	P04088
GeneID	397490;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: INHBB

Note: This product is for in vitro research use only