

Kinesin heavy chain isoform 5A

Catalog No: #AP84548



Package Size: #AP84548-1 50ug #AP84548-2 100ug #AP84548-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kinesin heavy chain isoform 5A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Kinesin heavy chain neuron-specific 1, Neuronal kinesin heavy chain, NKHC, NKHC1
Accession No.	Uniprot ID:Q12840
Uniprot	Q12840
GeneID	3798;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: KIF5A

Note: This product is for in vitro research use only