

# Interleukin-4 receptor subunit alpha

Catalog No: #AP84551



Package Size: #AP84551-1 50ug #AP84551-2 100ug #AP84551-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Interleukin-4 receptor subunit alpha                   |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Hu   |
| Other Names        | IL-4 receptor subunit alpha,CD124,IL4RA,582J2.1        |
| Accession No.      | Uniprot ID:P24394                                      |
| Uniprot            | P24394   |
| GeneID             | 3566;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL4R

Note: This product is for in vitro research use only