

Hypoxia-inducible factor 3-alpha

Catalog No: #AP84557



Package Size: #AP84557-1 50ug #AP84557-2 100ug #AP84557-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hypoxia-inducible factor 3-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HIF-3-alpha,HIF3-alpha-1
Accession No.	Uniprot ID:Q9JHS2
Uniprot	Q9JHS2
GeneID	64345;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hif3a

Note: This product is for in vitro research use only