

Major histocompatibility complex class I-related gene protein

Catalog No: #AP84568

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84568-1 50ug #AP84568-2 100ug #AP84568-3 1mg

Description

Product Name	Major histocompatibility complex class I-related gene protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MHC class I-related gene protein, Class I histocompatibility antigen-like protein, Mr1a
Accession No.	Uniprot ID:Q8HWB0
Uniprot	Q8HWB0
GeneID	15064;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mr1

Note: This product is for in vitro research use only