

## Interleukin-31 receptor subunit alpha

Catalog No: #AP84571



Package Size: #AP84571-1 50ug #AP84571-2 100ug #AP84571-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-31 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-31 receptor subunit alpha, GLM-R, mGLM-R, Gp130-like monocyte receptor, Gp130-like receptor, Novel cytokine receptor 10, NR10, ZcytoR17, Glmr
Accession No.	Uniprot ID:Q8K5B1
Uniprot	Q8K5B1
GeneID	218624;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: IL31ra

Note: This product is for in vitro research use only