

# Interleukin-36 beta

Catalog No: #AP84572



Package Size: #AP84572-1 50ug #AP84572-2 100ug #AP84572-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-36 beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Interleukin-1 family member 8,IL-1F8,Fil1e,Il1f8
Accession No.	Uniprot ID:Q9D6Z6
Uniprot	Q9D6Z6
GeneID	69677;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Il36b

Note: This product is for in vitro research use only