

## Histone acetyltransferase KAT8

Catalog No: #AP84574



Package Size: #AP84574-1 50ug #AP84574-2 100ug #AP84574-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Histone acetyltransferase KAT8   |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Ms   |
| Other Names        | Lysine acetyltransferase 8, MOZ, YBF2/SAS3, SAS2 and TIP60 protein 1, MYST-1, Mof, Myst1 |
| Accession No.      | Uniprot ID:Q9D1P2  |
| Uniprot            | Q9D1P2   |
| GeneID             | 67773;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                   |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Kat8

Note: This product is for in vitro research use only