

Inter-alpha-trypsin inhibitor heavy chain H4

Catalog No: #AP84579



Package Size: #AP84579-1 50ug #AP84579-2 100ug #AP84579-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inter-alpha-trypsin inhibitor heavy chain H4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	ITI heavy chain H4, Inter-alpha-trypsin inhibitor family heavy chain-related protein, IHRP, Major acute phase protein, MAP, IHRP
Accession No.	Uniprot ID:P79263
Uniprot	P79263
GeneID	396799;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ITIH4

Note: This product is for in vitro research use only