

Transcription factor AP-1

Catalog No: #AP84580



Package Size: #AP84580-1 50ug #AP84580-2 100ug #AP84580-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transcription factor AP-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Proto-oncogene c-Jun
Accession No.	Uniprot ID:P18870
Uniprot	P18870
GeneID	424673;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:JUN

Note: This product is for in vitro research use only