

Interleukin-13 receptor subunit alpha-2

Catalog No: #AP84584



Package Size: #AP84584-1 50ug #AP84584-2 100ug #AP84584-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-13 receptor subunit alpha-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-13 receptor subunit alpha-2,CD213a2
Accession No.	Uniprot ID:O88786
Uniprot	O88786
GeneID	16165;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL13ra2

Note: This product is for in vitro research use only