

E3 ubiquitin-protein ligase HACE1

Catalog No: #AP84589



Package Size: #AP84589-1 50ug #AP84589-2 100ug #AP84589-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase HACE1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HECT domain and ankyrin repeat-containing E3 ubiquitin-protein ligase 1,HECT-type E3 ubiquitin transferase HACE1
Accession No.	Uniprot ID:D3ZBM7
Uniprot	D3ZBM7
GeneID	361866;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hace1

Note: This product is for in vitro research use only