

C-Jun-amino-terminal kinase-interacting protein 4

Catalog No: #AP84592



Package Size: #AP84592-1 50ug #AP84592-2 100ug #AP84592-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-Jun-amino-terminal kinase-interacting protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	JIP-4, JNK-associated leucine-zipper protein, JLP, JNK/SAPK-associated protein 2, JSAP2, Mitogen-activated protein kinase 8-interacting protein 4, Sperm-associated antigen 9, Jip4, Jsap2, Kiaa0516, Mapk8ip4
Accession No.	Uniprot ID:Q58A65
Uniprot	Q58A65
GeneID	70834;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Spag9

Note: This product is for in vitro research use only