

Insulin receptor substrate 4

Catalog No: #AP84595



Package Size: #AP84595-1 50ug #AP84595-2 100ug #AP84595-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin receptor substrate 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IRS-4,160 kDa phosphotyrosine protein,py160,Phosphoprotein of 160 kDa,pp160
Accession No.	Uniprot ID:O14654
Uniprot	O14654
GeneID	8471;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IRS4

Note: This product is for in vitro research use only