

# Indolethylamine N-methyltransferase

Catalog No: #AP84598



Package Size: #AP84598-1 50ug #AP84598-2 100ug #AP84598-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Indolethylamine N-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Indolamine N-methyltransferase, Aromatic alkylamine N-methyltransferase, Amine N-methyltransferase, Thioether S-methyltransferase, TEMT
Accession No.	Uniprot ID:O95050
Uniprot	O95050
GeneID	11185;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: INMT

Note: This product is for in vitro research use only