

STMN2 (Phospho-Ser73) Antibody

Catalog No: #12644



Package Size: #12644-1 50ul #12644-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	STMN2 (Phospho-Ser73) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	STMN2 (Phospho-Ser73) Antibody detects endogenous levels of STMN2 only when phosphorylated at Ser73
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human STMN2 (Phospho-Ser73)
Target Name	STMN2
Modification	Phospho
Other Names	STMN2, SCGN10, SCG10, Stathmin-like 2, Protein SCG10, Stathmin-2
Accession No.	Swiss-Prot#: Q93045NCBI Gene ID: 11075
Uniprot	Q93045
GeneID	11075;
Target Species	human
Calculated MW	21kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only