

# Heat shock 70 kDa protein 1A/1B

Catalog No: #AP84602



Package Size: #AP84602-1 50ug #AP84602-2 100ug #AP84602-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Heat shock 70 kDa protein 1A/1B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Heat shock 70 kDa protein 1/2,HSP70-1/HSP70-2,HSP70.1/HSP70.2,Hsp70-1,Hspa1
Accession No.	Uniprot ID:Q07439
Uniprot	Q07439
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Hspa1a

Note: This product is for in vitro research use only