

## Heat shock cognate 71 kDa protein

Catalog No: #AP84603



Package Size: #AP84603-1 50ug #AP84603-2 100ug #AP84603-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Heat shock cognate 71 kDa protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Heat shock 70 kDa protein 8,Hsc70,Hsc73
Accession No.	Uniprot ID:P63017
Uniprot	P63017
GeneID	15481;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Hspa8

Note: This product is for in vitro research use only