

Keratin, type I cytoskeletal 13

Catalog No: #AP84604

Package Size: #AP84604-1 50ug #AP84604-2 100ug #AP84604-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Keratin, type I cytoskeletal 13 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Cytokeratin-13,CK-13,Keratin-13,K13,Type I keratin Ka13,Ka13 |
| Accession No. | Uniprot ID:Q61FV4 |
| Uniprot | Q61FV4 |
| GeneID | 287699; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Krt13

Note: This product is for in vitro research use only