

Keratin, type I cytoskeletal 9

Catalog No: #AP84605

Package Size: #AP84605-1 50ug #AP84605-2 100ug #AP84605-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytokeratin-9,CK-9,Keratin-9,K9
Accession No.	Uniprot ID:P35527
Uniprot	P35527
GeneID	3857;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT9

Note: This product is for in vitro research use only