

## Keratin, type II cytoskeletal 7

Catalog No: #AP84606



Package Size: #AP84606-1 50ug #AP84606-2 100ug #AP84606-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Keratin, type II cytoskeletal 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokeratin-7,CK-7,Keratin-7,K7,Type-II keratin Kb7,Kb7,Krt2-7
Accession No.	Uniprot ID:Q6IG12
Uniprot	Q6IG12
GeneID	300242;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Krt7

Note: This product is for in vitro research use only