

Insulin-like growth factor 1 receptor

Catalog No: #AP84610



Package Size: #AP84610-1 50ug #AP84610-2 100ug #AP84610-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like growth factor 1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Insulin-like growth factor I receptor,IGF-I receptor,CD221
Accession No.	Uniprot ID:P08069
Uniprot	P08069
GeneID	3480;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IGF1R

Note: This product is for in vitro research use only