

## Interleukin-22

Catalog No: #AP84613



Package Size: #AP84613-1 50ug #AP84613-2 100ug #AP84613-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-22
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-22,IL-10-related T-cell-derived-inducible factor,IL-TIF,IL-TIF alpha,Interleukin-22a,IL-22a,Il22a,Il22f,Il22g
Accession No.	Uniprot ID:Q9JJY9
Uniprot	Q9JJY9
GeneID	50929;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Il22

Note: This product is for in vitro research use only