

# Interleukin-4 receptor subunit alpha

Catalog No: #AP84616



Package Size: #AP84616-1 50ug #AP84616-2 100ug #AP84616-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-4 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	IL-4 receptor subunit alpha,CD124
Accession No.	Uniprot ID:Q863Z5
Uniprot	Q863Z5
GeneID	397614;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL4R

Note: This product is for in vitro research use only