

Interleukin-5 receptor subunit alpha

Catalog No: #AP84617



Package Size: #AP84617-1 50ug #AP84617-2 100ug #AP84617-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-5 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-5 receptor subunit alpha,CDw125,CD125,IL5R
Accession No.	Uniprot ID:Q01344
Uniprot	Q01344
GeneID	3568;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL5RA

Note: This product is for in vitro research use only