

Linker for activation of T-cells family member 1

Catalog No: #AP84620



Package Size: #AP84620-1 50ug #AP84620-2 100ug #AP84620-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Linker for activation of T-cells family member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	36 kDa phospho-tyrosine adapter protein,pp36,p36-38
Accession No.	Uniprot ID:O43561
Uniprot	O43561
GeneID	27040;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LAT

Note: This product is for in vitro research use only