

Serine/threonine-protein kinase LATS1

Catalog No: #AP84621



Package Size: #AP84621-1 50ug #AP84621-2 100ug #AP84621-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase LATS1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Large tumor suppressor homolog 1,WARTS protein kinase,h-warts,WARTS
Accession No.	Uniprot ID:O95835
Uniprot	O95835
GeneID	9113;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LATS1

Note: This product is for in vitro research use only