

6-phosphofructokinase, liver type

Catalog No: #AP84622



Package Size: #AP84622-1 50ug #AP84622-2 100ug #AP84622-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	6-phosphofructokinase, liver type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Phosphofructo-1-kinase isozyme B,PFK-B,Phosphofructokinase 1,Phosphohexokinase,Pfk-I,Pfkb
Accession No.	Uniprot ID:P12382
Uniprot	P12382
GeneID	18641;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PfkI

Note: This product is for in vitro research use only