

## Insulin-like growth factor II

Catalog No: #AP84623



Package Size: #AP84623-1 50ug #AP84623-2 100ug #AP84623-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Insulin-like growth factor II
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IGF-II, Multiplication-stimulating polypeptide, IGF-2
Accession No.	Uniprot ID:P09535
Uniprot	P09535
GeneID	16002;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: IGF2

Note: This product is for in vitro research use only