

DNA-binding protein Ikaros

Catalog No: #AP84624



Package Size: #AP84624-1 50ug #AP84624-2 100ug #AP84624-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA-binding protein Ikaros
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ikaros family zinc finger protein 1,Lymphoid transcription factor LyF-1,IK1,IKAROS,LYF1,ZNFN1A1
Accession No.	Uniprot ID:Q13422
Uniprot	Q13422
GeneID	10320;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IKZF1

Note: This product is for in vitro research use only