

Calcium/calmodulin-dependent protein kinase kinase

2

Catalog No: #AP84627

Package Size: #AP84627-1 50ug #AP84627-2 100ug #AP84627-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcium/calmodulin-dependent protein kinase kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CaM-KK 2, Calcium/calmodulin-dependent protein kinase kinase beta, CaM-KK beta, Kiaa0787
Accession No.	Uniprot ID:Q8C078
Uniprot	Q8C078
GeneID	207565;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Camkk2

Note: This product is for in vitro research use only