

Integrin alpha-9

Catalog No: #AP84630



Package Size: #AP84630-1 50ug #AP84630-2 100ug #AP84630-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Integrin alpha-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Integrin alpha-RLC
Accession No.	Uniprot ID:Q13797
Uniprot	Q13797
GeneID	3680;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ITGA9

Note: This product is for in vitro research use only