

Glandular kallikrein-7, submandibular/renal

Catalog No: #AP84634



Package Size: #AP84634-1 50ug #AP84634-2 100ug #AP84634-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glandular kallikrein-7, submandibular/renal
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	rGK-7, Esterase B, Kallikrein-related protein K1, Proteinase A, RSKG-7, Tissue kallikrein, Kik-7
Accession No.	Uniprot ID:P36373
Uniprot	P36373
GeneID	24523;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Kik7

Note: This product is for in vitro research use only