

Serine/threonine-protein kinase LATS2

Catalog No: #AP84637



Package Size: #AP84637-1 50ug #AP84637-2 100ug #AP84637-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase LATS2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Kinase phosphorylated during mitosis protein, Large tumor suppressor homolog 2, Serine/threonine-protein kinase kpm, Warts-like kinase, KPM
Accession No.	Uniprot ID:Q9NRM7
Uniprot	Q9NRM7
GeneID	26524;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LATS2

Note: This product is for in vitro research use only