

L-selectin

Catalog No: #AP84639



Package Size: #AP84639-1 50ug #AP84639-2 100ug #AP84639-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-selectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD62 antigen-like family member L,Leukocyte adhesion molecule 1,LAM-1,Leukocyte-endothelial cell adhesion molecule 1,LECAM1,Lymph node homing receptor,Lymphocyte antigen 22,Ly-22,Lymphocyte surface MEL-14 antigen,CD62L,Lnhr,Ly-22,Ly22
Accession No.	Uniprot ID:P18337
Uniprot	P18337
GeneID	20343;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sell

Note: This product is for in vitro research use only