

Protein-lysine 6-oxidase

Catalog No: #AP84640



Package Size: #AP84640-1 50ug #AP84640-2 100ug #AP84640-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein-lysine 6-oxidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Lysyl oxidase
Accession No.	Uniprot ID:P16636
Uniprot	P16636
GeneID	24914;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lox

Note: This product is for in vitro research use only