

Interferon-induced protein 44-like

Catalog No: #AP84641



Package Size: #AP84641-1 50ug #AP84641-2 100ug #AP84641-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon-induced protein 44-like
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	H28
Accession No.	Uniprot ID:Q9BDB7
Uniprot	Q9BDB7
GeneID	15061;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ifi44l

Note: This product is for in vitro research use only