

L-lactate dehydrogenase C chain

Catalog No: #AP84647



Package Size: #AP84647-1 50ug #AP84647-2 100ug #AP84647-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-lactate dehydrogenase C chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LDH-C,Cancer/testis antigen 32,CT32,LDH testis subunit,LDH-X,LDH3,LDHX
Accession No.	Uniprot ID:P07864
Uniprot	P07864
GeneID	3948;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LDHC

Note: This product is for in vitro research use only