

# Transcription factor AP-1

Catalog No: #AP84649



Package Size: #AP84649-1 50ug #AP84649-2 100ug #AP84649-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transcription factor AP-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Activator protein 1,AP1,Proto-oncogene c-Jun,V-jun avian sarcoma virus 17 oncogene homolog,Rjg-9
Accession No.	Uniprot ID:P17325
Uniprot	P17325
GeneID	24516;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Jun

Note: This product is for in vitro research use only