

Laminin subunit beta-2

Catalog No: #AP84658



Package Size: #AP84658-1 50ug #AP84658-2 100ug #AP84658-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Laminin subunit beta-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Laminin chain B3,Laminin-11 subunit beta,Laminin-14 subunit beta,Laminin-15 subunit beta,Laminin-3 subunit beta,Laminin-4 subunit beta,Laminin-7 subunit beta,Laminin-9 subunit beta,S-laminin subunit beta,S-LAM beta
Accession No.	Uniprot ID:P15800
Uniprot	P15800
GeneID	25473;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lamb2

Note: This product is for in vitro research use only