

## Insulin-like growth factor 1 receptor

Catalog No: #AP84662



Package Size: #AP84662-1 50ug #AP84662-2 100ug #AP84662-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Insulin-like growth factor 1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Insulin-like growth factor I receptor, IGF-I receptor, CD221
Accession No.	Uniprot ID:Q05688
Uniprot	Q05688
GeneID	281848;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: IGF1R

Note: This product is for in vitro research use only