

Interleukin-1 alpha

Catalog No: #AP84665



Package Size: #AP84665-1 50ug #AP84665-2 100ug #AP84665-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	IL-1 alpha
Accession No.	Uniprot ID:O46612
Uniprot	O46612
GeneID	403782;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL1A

Note: This product is for in vitro research use only