

Interleukin-1 receptor-associated kinase 1

Catalog No: #AP84666



Package Size: #AP84666-1 50ug #AP84666-2 100ug #AP84666-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 receptor-associated kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IRAK,Pelle-like protein kinase,mPLK,Il1rak
Accession No.	Uniprot ID:Q62406
Uniprot	Q62406
GeneID	16179;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Irak1

Note: This product is for in vitro research use only